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**EN Wellion® CALLA test strips** IMPORTANT: PLEASE READ THIS INFORMATION AND YOUR WELLION CALLA USER GUIDE BEFORE USING THE CALLA TEST STRIPS TO TEST YOUR BLOOD SUGAR.

Please call our authorized representative by using the information printed at the end of this insert or on the box or visit our website.

**INTENDED USE** The Wellion CALLA Blood Glucose Test Strips are to be used with the Wellion CALLA Blood Glucose Meters to quantitatively measure glucose in capillary whole blood. The Wellion CALLA System is plasma-calibrated for easy comparison to lab result. The Wellion CALLA Blood Glucose Monitoring System is intended for self-testing by persons with diabetes and by health care professionals.

**Introduction** The biosensor-based Wellion CALLA Blood Glucose Test Strips employs the enzyme Glucose Oxidase that is specific for glucose measurement from 20 to 600 mg/dL (1,1 to 33,3 mmol/L). The arrow printed on the test strip indicates the side of the test strip where it slips in the blood sample. When blood from a puncture site touches the tip of the test strip, capillary action leads the blood into the reaction zone. The enzyme reacts with glucose in blood and produces an electrical signal. This signal is measured by the Wellion CALLA Meter. Test result will be displayed in 6 seconds.

**Reagent Composition** Each cm<sup>2</sup> of test strip contains the following ingredients in the approximate concentrations listed below:

Glucose Oxidase (A. Niger)	2,3%
Electron Shuttle	2,0%
Non-Reactive Ingredients	72,7%

**Warnings and Precautions** [1] The Wellion CALLA Blood Glucose Test Strips are for use outside the body (IN VITRO diagnostic use). [2] Do not use test strips after expiration date. The expiration date is printed on the carton and on the label of the test strip bottle. [3] Do not use strips that are torn or damaged. [4] Do not re-use the strips. [5] If your test result is below 50 mg/dL (2,8 mmol/L) or above 250 mg/dL (13,9 mmol/L), perform the Control Solution Test to make sure your system is working properly and then repeat testing once more using blood from your fingertip. If the result is still the same; contact your healthcare professional immediately. [6] If you have symptoms or are inconsistent with your test results, and you have eliminated common errors as described in the user's guide, contact your healthcare professional immediately. [7] Never make significant changes in your diabetes treatment program or ignore symptoms without consulting your physician. [8] Use the Wellion CALLA Blood Glucose Monitoring System within temperature range of 10°C to 40°C (50°F to 104°F).

**Storage and Handling** [1] Always close the vial cap tightly after removing a test strip to avoid moisture and direct sunlight. [2] Store the test strips at temperatures between 4°C to 30°C (39°F to 86°F). [3] Store out of direct sunlight; do not refrigerate or freeze. [4] Unopened test strips are stable until the expiration date printed on the box when stored properly. [5] The month printed refers to the END of that month. [6] Use within 6 months after first opening. [7] Do not handle the test strips with wet or dirty hands.

**Sample Collection and Preparation** Fresh finger stick samples that blood may be collected into heparin or sodium EDTA test tubes and should be used within 15 minutes. Do not test on plasma or serum samples. Do not use fluoride-based preservatives.

**Test Procedure** Read "Test Your Blood Glucose" section in the User's Guide of the Wellion CALLA Blood Glucose Meter.

**Alternate Site Testing** The Wellion CALLA Blood Glucose Monitoring System can also test for blood glucose from palm and forearm (alternate site testing, or AST). Physiological differences between the finger or palm and forearm may result in different blood glucose readings. [1] You should not test on palm and forearm and use finger tip testing instead when:

- [1] You have hypoglycemic unawareness. [2] You suspect your blood glucose level is changing rapidly because of exercise or after rapid-acting insulin injection or it has been less than two hours after a meal. [3] You will be driving a car or operating any machinery activities such as operating machinery. [4] Consult your healthcare professional to decide if alternate site testing is right for you.

**Expected Values** Blood glucose value depends on many factors. The ideal range for adult with diabetes is:

- Less than 100 mg/dL (5,6 mmol/L) before meals
- Less than 140 mg/dL (7,8 mmol/L) after meals

It is important to consult with your physician or healthcare professional to determine an appropriate target range for you.

**Quality Control** Run Level 0, Level 1 or Level 2 control solution test. Follow the User's Guide instructions. You can perform control solution testing when:

- [1] Your test results do not agree with how you feel. [2] Once per week to verify that the meter and test strips are working properly together. [3] Your test strips were stored at temperatures and humidity outside the specified storage conditions. [4] When you open a new bottle of test strips. [5] Practice your testing technique. [6] You drop your bottle of test strips.

Use only Wellion CALLA Control Solution. Your test results should fall within the control range printed on the test strip bottle. If test results falls outside the control range, repeat the test. Results that fall outside the range may be caused by:

- [1] Error in performing the control solution test. [2] Expired or contaminated control solution. [3] Test Strip is damaged. [4] Meter malfunction.

If the result continues to fall out side of the printed range; please call our authorized representative by using the information printed at the end of this insert or on the box for assistance.

**Problem Solving** [1] Confirm the Test Strips are within the expiration date. [2] Make sure the blood completely fills the reaction zone. DO NOT add a second drop of blood; discard the test strip from the meter and using a new Test Strip. [3] Check the system by performing the Control Solutions Test. [4] Refer to the "Solving Problems" section in the User Guide for additional information. Please call our authorized representative by using the information printed at the end of this insert or on the box for assistance or your healthcare professional.

**Limitations** [1] DO NOT use plasma or serum samples. [2] DO NOT test on neonatal samples (new born). [3] DO NOT test on venous or arterial samples. [4] Altitude up to 2,000 meters will not affect test results. [5] Hematocrit range: 35% to 50%. [6] Triglycerides: 260 to 360 mg/dL (2,9 to 4,0 mmol/L) has no significant effect on test results (normal range 36 to 165 mg/dL or 2 to 9,2 mmol/L). [7] Icodextrin and its metabolites (maltoze, maltotriose and maltotetraose) do not significantly affect test results. [8] The following will not affect test results in expected blood concentrations: ascorbic acid (Vitamin C), uric acid, or methyl-dopa. [9] Therapeutic concentration of dopamine or L-dopa may affect results. Individuals taking these drugs should not use this system to test for blood sugar. [10] Higher than therapeutic concentrations of acetaminophen may affect test results. Individuals taking these drugs should not use this system to test for blood sugar.

**PERFORMANCE EVALUATION DATA** **Accuracy** According to EN ISO 15197:2015, total of 600 fingertip blood data were collected for system accuracy evaluation, 210 data were from glucose level below 100 mg/dL (5,6 mmol/L) and 390 data were from glucose level above or equal to 100 mg/dL (5,6 mmol/L). The system accuracy evaluation of fingertip were tested by healthcare professional, and compared to YSI 2300 Blood Glucose Analyzer.

Table 1: Accuracy for blood glucose level <100 mg/dL (5,6 mmol/L)

Izvor uzorka	Unutar ±5mg/dL (±0,3mmol/L)	Unutar ±10mg/dL (±0,6mmol/L)	Unutar ±15mg/dL (±0,8mmol/L)
Kapilarna krv	108/210 (51,4%)	167/210 (79,5%)	205/210(97,6%)

Table 2: Accuracy for blood glucose level ≥100 mg/dL (5,6 mmol/L)

Sample Source	Within ± 5%	Within ± 10%	Within ± 15%
Capillary	272/390 (69,7%)	366/390 (93,8%)	387/390 (99,2%)

Table 3: Accuracy for blood glucose level <100 mg/dL (5,6 mmol/L)

Sample Source	Within ± 5mg/dL (±0,3mmol/L)	Within ± 10mg/dL (±0,6mmol/L)	Within ± 15mg/dL (±0,8mmol/L)
Capillary	116/210 (55,2%)	186/210 (88,6%)	210/210 (100%)

Table 4: Accuracy for blood glucose level ≥100 mg/dL (5,6 mmol/L)

Sample Source	Within ± 5% (Unutar ±10%)	Within ± 10% (Unutar ±15%)	
Capillary	262/390 (67,2%)	367/390 (94,1%)	386/390 (99,0%)

Table 5: Accuracy for blood glucose level <100 mg/dL (5,6 mmol/L)

Sample Source	Within ± 15mg/dL (±1,5%)
Capillary	596/600 99,3%

Table 6: Accuracy for blood glucose level <100 mg/dL (5,6 mmol/L)

Sample Source	Within ± 5mg/dL (±0,3mmol/L)	Within ± 10mg/dL (±0,6mmol/L)	Within ± 15mg/dL (±0,8mmol/L)
Capillary	110/210 (52,4%)	175/210 (83,3%)	207/210 (98,6%)

Table 7: Accuracy for blood glucose level ≥100 mg/dL (5,6 mmol/L)

Sample Source	Within ± 5% (Unutar ±10%)	Within ± 10% (Unutar ±15%)	
Capillary	252/390 (64,6%)	353/390 (90,5%)	387/390 (99,2%)

Table 8: Accuracy for blood glucose level <100 mg/dL (5,6 mmol/L)

Sample Source	Within ± 15mg/dL (±1,5%)
Capillary	596/600 99,2%

The result showed the Wellion CALLA System compares well with the YSI 2300 Glucose Analyzer and meets the minimum acceptable performance criteria defined in EN ISO 15197:2015.

**Precision** Three lots of Wellion CALLA Test Strips were used for within-run precision study. Venous blood in heparin-tubes was spiked to 5 concentrations. Blood glucose readings were recorded for 1 day resulting in 300 data points for each concentration; as shown in the following tables:

With Wellion CALLA Dialog Meter

Number of Readings	300	300	300	300
Average Reading mg/dL (mmol/L)	50 (2,8)	71 (3,9)	134 (2,2)	